Specifications

The state of the s	REFRACTIVE POWER MEASUREM	
	Spherical refractive power (S)	-30.00 D to +25.00 l (at VD = 12.0 mm)
	Cylindrical power (C)	0.00 D to ±12.50 D (at VD = 12.0 mm)
	Astigmatic axis (A)	0° to 180°
	KERATOMETRY MEASUREMENT	
	Corneal curvature radius	5.00 mm to 13.00 m
	Corneal astigmatic axis	0° to 180°
	INTRAOCULARRESSUREMEASUR	EMENT
	Measurement range	1 mmHg to 60 mmH (1 hPa to 80 hPa)
	PACHYMETRY MEASUREMENT	
	Measurement range	300 µm to 800 µm
	TOPOGRAPHY MEASUREMENT	
	Corneal curvature radius	5.50 mm to 10.00 m
	Corneal astigmatic axis	0° to 180°
	AUXILIAR¥UNCTION	
	Interpupillary distance	Measurement range 20 mm to 85 mm
	Corneal diameter and pupil diameter	Measurement range 1.00 mm to 14.00 m
	Dry Eye observation	Blinking analysis, Tear meniscus heigł Hyperaemia, Meibomian glands

Interna

Printe

Data o


Export


al database	Integrated SD card
	Integrated thermal printer
utput type	3 × USB 2.0, 1 x LAN, 1 × SD card slot, 1 × WLAN
format	CSV, XML, JPG, PDF

Dimensions WDH	$312 \times 491 \times 450 \text{ mm}$
Veight	Approx. 23 kg
/oltage	100 VAC to 240 VAC
Frequency	50/60 Hz
Power consumption	110 VA

Dry Eye observation software

The combination of 4 examinations gives a detailed impression of the patients' ocular surface. Observe the blinking frequency, tear meniscus height, Hyperemia and Meibomian glands by using only one medical device.





Blinking frequency

Tear meniscus height







Hyperemia


Meibomian glands

Topography analysis software





Dual map instantaneous





Single map

Fourier analysis





RODENSTOCK Instruments



The all-new **ALINO[®] MULTIFUNCTION UNIT**

Making analysis fashionable



The all-new ALINO Tune up your business

Experience 6 analyses in 1 device

Are you ready for something that suits your business perfectly? Let's have a look at how powerful ALINO[®] can be.

What's your benefit?

Be competitive and provide added value compared to online stores and competitors.



6 in 1 combination Impress your custom ers with this ALL-IN-ONE device

More efficient

11

Higher customer

circulation = increased

frequency of customers



Time saving



Auto alignment and auto measurement

Generate regular additional saleand retain customers for the long-term

Cost reduction The combination of all product features reduces your investment while providing added value



customer journey! Thanks to the high amount of examinations and analyses, you gain a better informative value for the client. This leads to increased consulting competence and stronger customer loyalty. It also creates trust and your customers will be more satisfied.



The ALINOmakes my daily work run much smoother! More benefit for my consumers & less effort for me!



ALIN^Osimplifies your workflow.



6 in 1 solution

01 Refraction

Having a good starting value for subjective refraction is essential. The QUICK REF MODE supports this even for uncooperative patients.

02 Keratometry

KAI (Kerato-asymmetry index) & KRI (Kerato-regularity index) give the first hint of possible irregularities on the cornea. This supports evaluation of the best vision correction or best-fitting contact lens.

03 Tonometry

The very gentle air puff created by a new generation of air-flow technology means that patients more readily accept the need for regular check-ups of the **IOP** measurement.

06 Dry Eye

Observe the customers' eye health and evaluate the long-term tolerability of contact lenses.

05 Topography

Various topographic maps support you in carrying out vision screenings, contact lense fittings and in patient education.

04 Pachymetry

Correcting the IOP by measuring the central corneal thickness is essential in interpreting the intraocular pressure. Take advantage of the combination of tests:

Additional services Tonometry + Pachymetry

Contact lens fitting and tolerance control Keratometry + Topography + Dry Eye

Eve healthcare Topography + Tonometry + Pachymetry + Dry Eye



rodenstock-instruments.de